







The LFGI Series provides either positive or negative pressure HEPA-filtered unidirectional (laminar) airflow for patient, personnel, and product protection. The unit can be configured for positive (non-hazardous) or negative (hazardous) pressure to comply with USP <797> and NIOSH recommendations.

<u> CERNIERE</u>

GERMFREE.com

800.888.5357


GERMFREE.com 800.888.5357

All Stainless Steel Laminar Flow Glovebox Isolator

Room Air

Filtered Air

- Meets USP <797> compliance.
- All surfaces, both inside and outside, are constructed from stainless steel to provide exceptional durability and strength.
- Unidirectional (laminar) airflow in the work area maintains ISO Class 5/Class 100 air quality during manipulations.
- Sliding work trays facilitate easy cleaning and simplify movement within the Isolator reducing operator reach and strain.
- Antechamber can be installed on either left, right, or both sides of the unit.
- Ergonomic placement of sharps and waste tubes.
- Programmable, height adjustable, ergonomic stand.
- Gloveports have a 3" armrest to reduce operator fatigue.
- Prefilter extends the life of the HEPA filters and can be replaced without tools.
- All corners in work area, and work surfaces are easily reached and cleaned.
- UL Listing 61010-1.
- Designed to fit through standard door openings and elevators.
- Extra large oval gloveports are placed with bottoms together to provide an anthropometrically correct configuration that accommodates a wide range of body types and increases range of movement.



configuration



Positive configuration

Equipment Specifications:

Overall Dimensions:

LFGI-3USP: 54"W (39"W at base) x 32"D x 80"H LFGI-4USP: 66"W (50"W at base) x 32"D x 80"H LFGI-6USP: 90"W (74"W at base) x 32"D x 80"H

Sealed two-door airlock:

(16"W x 13 1/2"D x 14 1/4"H)

Work Area Dimensions:

LFGI-3USP: 34 1/2" W x 24 1/2" D x 29" H LFGI-4USP: 46 1/2" W x 24 1/2" D x 29" H LFGI-6USP: 70 1/2" W x 24 1/2" D x 29" H





Performance Data:

The Laminar Flow Glovebox/Isolator (LFGI) provides the ultimate in both product and operator protection. The LFGI functions as a Glovebox while operating under negative pressure to meet NIOSH recommendations and as an Isolator under positive pressure to meet USP <797> regulations. HEPA- filtered unidirectional (laminar) airflow bathes the work area to protect the product from contamination and removes any generated particulates. The operator is provided a complete barrier from materials being manipulated in the work area. The LFGI Series of barrier isolators meet or exceed ISO Class 5/Class 100 air quality. Each LFGI, before shipping, undergoes rigorous physical testing to assure the unit meets performance requirements. It is required that independent certification be performed before use.

General Information:

Work trays slide over one another and are removable to facilitate easy clean-up under the work deck. The large acrylic viewing panel is encased in a stainless steel frame and is at an angle to reduce glare and operator strain. Also, the viewing panel is top hinged and self supporting for easy access to the work area for cleaning and equipment ingress. Straight sides and back maximize work area to accommodate the many types and shapes of pharmacy equipment and compounders. Manual or programmable electric hydraulic assist height-adjusting stand offers a full 10" range variance allowing operators to sit or stand comfortably for extended time periods.

Antechamber:

Airlock maintains complete environmental separation between work area and ambient conditions. HEPA-filtered purge of trapped airlock air eliminates cross-contamination between the work area and the room during both material ingress and egress.

Waste Removal:

Trash tube is a 4" diameter tube, flanged for waste bag/container connection and provided with an easily removable seal/stopper to maintain work area containment. Sharps tube is a straight 2 1/2" diameter tube to facilitate quick drop of the largest syringes and is provided with an easily removable seal/stopper to maintain work area containment.

Configuration Options:

Overall work area pressure can be either positive or negative and can be altered in place by a certifier. Antechamber can be left or right handed or can be installed on both sides. Waste and sharps removal tubes can be left or right handed or can be installed on both sides.

Filtration:

Full framed minipleat HEPA- filters 100% of the inlet and exhaust air supplied to the work area and the antechamber. HEPA filters are front loading for easy replacement by a certifier.

Exhaust:

Clean air plenum is external to the work area. Constructed to allow for optional outside venting of exhaust air. Venting is not required for operation but is optional.



- Complete the contact information
- 2 Check the boxes next to the desired products
- **3** Fax this form to 386-677-1114

Name: _

4 A GERMFREE specialist will contact you

If you have questions, please call 800-888-5357.

	Company:		
	Address:		
	City/State/Zip:		
	Phone:()	Ext:Email:	
	Fax:()		
	GPO contract #:		
All products available under contract with GSA.			
Qty.	Part Number	Description	
()	3 Foot Laminar Flov LFGI-3USP	w Glovebox Isolator 3 Foot Laminar Flow Glovebox Isolator and Accessories	
()	4 Foot Laminar Flow LFGI-4USP	w Glovebox Isolator 4 Foot Laminar Flow Glovebox Isolator and Accessories	
()	6 Foot Laminar Flow Glovebox Isolator LFGI-6USP 6 Foot Laminar Flow Glovebox Isolator and Accessories		
	Configuration Choices for all LFGI Series Units Overall Pressure: Airlock placement: Right Sharps/Waste placement: Right Right Left		
	Optional Accessories for all LFGI Series Units		
	LFGI-A12 LFGI-A33B LFGI-A33TGOS LFGI-A24 LFGI-A402C LFGI-A51D LFGI-A80B LFGI-A83 LFGI-A85M LFGI-A85S LFGI-A8E GRX	Exhaust Transition Additional Pair of Smooth Butadyl/Hypalon Sleeves	